



SORATH ORTHO

A DIVISION OF SIGMA SURGICAL

HYPER INTERLOCKING NAIL FOR TIBIAL BONE

Providing Solutions in Orthopaedics...



ISO 13485

www.sorathortho.co.in


Hyper Interlocking Nail For Tibial Bone

Hyper Interlocking Nail For Tibial Bone

Dia	Size	Titanium	Stainless Steel
8mm	280mm to 380mm	E01304280 to E01304380	E01301280 to E01301380
9mm	280mm to 380mm	E01404280 to E01404380	E01401280 to E01401380
10mm	280mm to 380mm	E01504280 to E01504380	E01501280 to E01501380
11mm	280mm to 380mm	E01604280 to E01604380	E01601280 to E01601380

(20mm Variation in Long Nail Sizes)

Note : 400mm Nail in all dia manufacture against order only.



Interlocking Bone Screw, Dia 3.9 mm

Size	Titanium	Stainless steel
26mm to 60mm	A04904026 to A04904060	A04901026 to A04901060

*(use for 8mm Nail Distal Hole) (26mm to 50mm, 2mm Variation in size)
(55mm to 60mm, 5mm Variation in size)



Interlocking Bone Screw, Dia 4.9 mm

Size	Titanium	Stainless steel
28mm to 80mm	A05004028 to A05004080	A05001028 to A05001080


(28 mm to 50 mm, 2mm Variation in size)
(55 mm to 80 mm, 5mm Variation in size)



End Cap For Hyper Interlocking Nail For Tibial Bone

Size	Titanium	Stainless steel
00mm to 15mm	E13004000 to E13004015	E13001000 to E13001015

(5mm Variation in sizes)



Code No. J0200006 Hyper Interlocking Nail For Tibial Bone Complete Instrument Set



Sr. No	Code	Description
1	J0010001	Bone AWL with T Handle
2	J0020003	Conical Bolt
3	J0030001	Depth Guage
4	J0040003	Distal Device
5	J2200300	Drill Bit, 4.0 mm x 300 mm
6	J2180250	Drill Bit, 3.2 mm x 250 mm
7	J0060001	Drill Sleeve, 8.0 mm x 4.0 mm
8	J0060002	Drill Sleeve, 8.0 mm x 3.2 mm
9	J0080003	Guide Wire Pusher
10	J1900115	Guide Wire, 2.5 mm x 115 cm (simple)
11	J2010115	Guide Wire, 2.5 mm x 115 cm (Bullet Tips)
12	J0090002	Combined Hammer
13	J0100002	Hammering Bolt
14	J0110002	Impactor Extractor
15	J0120001	Protection Sleeve, 10.0 mm x 8.0 mm (Short)
16	J0130003	Main Zig
17	J0320001	Oblique Proximal Locking Device
18	J0140001	Screw Driver, 3.5 mm, Hexagonal (For 4.9mm interlocking Screw)
19	J0390001	Allen Key
20	J0160001	Teflon Sleeve
21	J0170001	Tissue Protector
22	J0180001	Trocar
23	J0570001	Ruler Scale
24	J0220001	Implant Tray
25	J1810002	Flexible Reamer Shaft
26	J1800002	Flexible Reamer Tip 8.5mm
27	J1800003	Flexible Reamer Tip 9mm
28	J1800004	Flexible Reamer Tip 9.5mm
29	J1800005	Flexible Reamer Tip 10mm
30	J1800006	Flexible Reamer Tip 10.5mm
31	J1800007	Flexible Reamer Tip 11mm
32	J1800008	Flexible Reamer Tip 11.5mm
33	J0190001	Instrument Box
34	J0200006	Complete Instrument Set For Hyper Tibial Nail

end to end

GLOBAL MEDICAL SUPPLY SOLUTIONS

○ Sourcing ○ Production ○ Quality Control ○ Logistics ○ Customer Services



- * Hyper, Sorath, Lexa, Spyrox, Sigma, Pelviflex and Flexora are trademarks or registered trademarks of SIGMA SURGICAL PVT. LTD. in INDIA, the USA and other Countries.
- * The information in the Catalogue is valid as of OCTOBER 2024 Designs and specifications are subject to be changed without notice.
- * Product shown in Catalogue may differ from actual products.
- * Sigma Surgical is not responsible for differences between the information in the Catalogue and the actual implants.

Manufacturing Facility :

SIGMA SURGICAL PRIVATE LIMITED

Survey no:1636/1,Opp. Bharat petrol pump,
Indore highway, Village-Kuha,
Taluka-Daskroi, Ahmedabad-382433.
Gujarat, India. 📍: 22.99946, 72.78633

Customer care :

+91-7575011006 / +91-8866211414

www.sorathortho.com
info@sorathortho.com