




SORATH ORTHO

A DIVISION OF SIGMA SURGICAL

Twin Lock-BCP Distal Femoral Plate

Twin Lock - Bone Compression Plate System



ISO 13485

www.sorathortho.co.in

FEMUR DISTAL

Twin Lock-BCP Distal Femoral Plate



Primary use : Distal Femur Fixation
Screw Size : Twin Lock-BCP Screw \varnothing 5.0mm/ Twin Lock-BCP Cancellous Screw \varnothing 5.0mm / Cortical Bone Screw \varnothing 4.5mm/ Twin Lock-BCP Cancellous Screw, \varnothing 5.0 mm
Matching Instrument : J0200020 Complete Instrument Set for 5.0 mm/ 4.5 Mm Locking System

Size (Holes)	Reference Number		Size (Holes)	Reference Number	
	Titanium	Stainless Steel		Titanium	Stainless Steel
	Left			Right	
5	C00314005	C00311005	5	C00324005	C00321005
7	C00314007	C00311007	7	C00324007	C00321007
9	C00314009	C00311009	9	C00324009	C00321009
11	C00314011	C00311011	11	C00324011	C00321011
13	C00314013	C00311013	13	C00324013	C00321013
15	C00314015	C00311015	15	C00324015	C00321015

end to end

GLOBAL MEDICAL SUPPLY SOLUTIONS

○ Sourcing ○ Production ○ Quality Control ○ Logistics ○ Customer Services



- * Hyper, Sorath, Lexa, Spyrox, Sigma, Pelviflex and Flexora are trademarks or registered trademarks of SIGMA SURGICAL PVT. LTD. in INDIA, the USA and other Countries.
- * The information in the Catalogue is valid as of OCTOBER 2024 Designs and specifications are subject to be changed without notice.
- * Product shown in Catalogue may differ from actual products.
- * Sigma Surgical is not responsible for differences between the information in the Catalogue and the actual implants.

Manufacturing Facility :

SIGMA SURGICAL PRIVATE LIMITED

Survey no:1636/1,Opp. Bharat petrol pump,
Indore highway, Village-Kuha,
Taluka-Daskroi, Ahmedabad-382433.
Gujarat, India. 📍 : 22.99946, 72.78633

Customer care :

+91-7575011006 / +91-8866211414

www.sorathortho.com
info@sorathortho.com